

## Folding Bed



The HT bed unit is integrated flush into the wall panel as a rest area in intensive care wards and patients' rooms. Since it can be quickly folded down and up, it is the ideal solution for relatives who wish to remain overnight. In an unused condition the

folding bed closes flush with the wall panels. In a used condition the folding bed provided an almost equivalent replacement for a normal bed.

### Dimensions

Width	2100 mm
Installation depth	252 mm
Height	1000 mm (in closed position)
Installation height	300 mm
Mattress	2000 x 900 x 80 mm

### Technical data

Bed frame	square profiles 80 x 40 x 2 mm; 40 x 40 x 2 mm; 20 x 20 x 2 mm made of polished stainless steel (Mat. no. 1.4301); welded together
Bed housing	rear and side panels made of 1.5 mm galvanised steel
Front panel	1 mm galvanised steel; with plaster board glued to the rear
Opening/closing	open position 90°: space in the room approx. 2090 x 800 mm closed position 0°: bed is flush with the wall panels
Handle	stainless steel bar handle
Bed fittings	2 x gas pressure string with 80 kg
Media rail (optional)	with electrical sockets for power and IT
Lighting unit (optional)	LED bed lamp

